

Original GE LogiqI LogiqE Transducer Linear Array Ultrasound Probe 12L-RS

Our Product Introduction

Basic Information

- Place of Origin: USA
- Brand Name: GE
- Model Number: 12L-RS
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 30-50pcs/month



Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Before Delivery: Tested Images&videos Will Be Shared
- Shipping Method: Express Or As Client's Required
- Highlight: **GE LogiqI Linear Array Ultrasound Probe, Linear Transducer Probe 12L-RS, Musculoskeletal Linear Transducer Probe**



More Images



for more products please visit us on 3dultrasoundprobe.com

Product Description

GE LogiqE LogiqE linear array original ultrasound transducer probe 12L-RS

1. Model: 12L-RS
2. Description: Peripheral Vascular, Small Parts, Musculoskeletal, Nerve Blocks, Superficial Thoracic/Pleural
3. Type: Linear array
4. Frequency: 5.0-13.0 Mhz
5. Compatible machine: Logiq e, Logiq i, Logiqbook, Logiqbook XP, Vivid i, Vivid q, Vivid S5, Vivid S6, Voluson e, Voluson i, Voluson P6, Voluson S6, Voluson S8





VOLKA MEDICAL TECH CO., LIMITED

+8619536805795

info@volka-medtech.com

3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong