


GE VividE Original Linear Ultrasound Probe 9L-RS , GE Ultrasound Transducers 3.3-10.0 Mhz

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com

Basic Information

- Brand Name: GE
- Model Number: 9L-RS
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 50pcs/month




Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Before Delivery: Tested Images&videos Will Be Shared
- Highlight: **3.3-10.0 Mhz Ge Ultrasound Transducers, GE VividE Linear Ultrasound Probe**



More Images




Product Description

GE linear array original ultrasound transducer probe 9L-RS

Type: Linear array
 Frequency: 3.3-10.0 Mhz
 Compatible machine: Logiq E, VIVID E
 Application: Peripheral Vascular, Vascular Access, Musculoskeletal
 Warranty: 60 days
 Condition: Used original one, fully tested before delivery

Note: Other probe models for GE Logiq/Vivid/Voluson series are available.



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong