

GE VolusonE6 VolusonE8 Ultrasound Repair Parts RFM423 Board KTZ304074

Our Product Introduction


Basic Information

- Place of Origin: USA
- Brand Name: GE
- Model Number: KTZ304074
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 working days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 60 Days After Fixed
- Service Time: 3-5 Working Days(normal Fault)
- Service: Repair
- Shipping Method: Express,or As Clients Required
- Highlight: **Ultrasound Accessories KTZ304074,
Ultrasound Repair Parts RFM423 Board,
GE Ultrasound Machine Spare Parts**



for more products please visit us on 3dultrasoundprobe.com

Product Description

GE VolusonE6 VolusonE8 Ultrasound Repair Parts RFM423 Board KTZ304074

Product Specifications

Attribute	Value
Warranty	60 days after fixed
Service time	3-5 working days (normal fault)
Service	Repair
Shipping method	Express, or as clients required

Product Description

Repair service for GE VolusonE6/VolusonE8 original RFM423 ultrasound board KTZ304074.

Part number: KTZ304074

Type: RFM423

Trade term: Repair service provided

Warranty: 60 days after fixed

Compatible Models

Repairs available for Voluson E6/E8/E10 (from BT15 to BT21)

RFM423 FE Mainboard Details

Front-End board for processing, generating, receiving and transmitting ultrasonic signals (used in Voluson E10 starting from BT17 to BT21).

Common Problems Solved

Hardware error DMA 0x408D (DMA buffer error)

Intense white interference on screen when sensor is turned on

Power button blinking yellow when turned on

Loading window appears without ultrasonic part (simulation mode activated)



More details welcome to contact with us!

Additional Parts in Stock

Brand	System	Description
GE	Logiq 200	MSTE board 2266030
GE	Logiq 200	mother board 2220786-2
GE	Logiq 400	PCONT(high voltage) 2269737
GE	Logiq 400	TLMC 2158390-2
GE	Logiq 400	VIPP 2245652
GE	Logiq 400	MPU ASSY (CPUboard) 2123339
GE	Logiq 400	LV3 power board 2211491
GE	Logiq 400	MASC board 2123317
GE	Logiq 400	TI board 2123324
GE	Logiq book	AC adapter 2396918
GE	Logiq book	MST mainboard 2365715

GE	Logiq book	DCDC ASSY 2382377
----	------------	-------------------

Frequently Asked Questions

Q: What is the MOQ?

A: 1pc.

Q: How's the stock condition?

A: Wide variety of Ultrasound probe/board, new&used&compatible new for your choose.

Q: Is there any guarantee for the goods?

A: Fully-tested and safe packing before shipping, less than 3% warranty return claims of previous orders.

Q: Will you provide warranty?

A: Items are sold with NO WARRANTY, unless it is stated otherwise in the body of the proforma invoice/purchase order. Mostly 90 days for new original, 30/60 days for used.

Q: Do you test all your goods before delivery?

A: Yes, we have 100% tested our goods before delivery.

Q: Can I have the goods returned if it is not working good?

A: Yes, goods are permitted to be returned in same received condition within 7 days. Refund will be handled within 2-4 working days after final confirmation (Freight & bank fee excluded). Shipping fee for return should be borne by buyers.

Q: If faults occurred after warranty date, can I contact you?

A: Sure, feel free to contact even years later. We'll look into the faults and try to find out some solutions/suggestions for you reference.

Q: How about the payment terms?

A: All orders are 100% prepaid before delivery. We accept wire transfer/ Western Union from bank, and cash. (US dollars, Euro, or CNY).

Q: How about delivery?

A: We ship worldwide. DHL/EMS/FedEx Express will be arranged within 1-5 working days after payment confirmed unless otherwise specified. We Seller are responsible for export Customs clearance in our end, and buyers for import Customs clearance in their end.



VOLKA MEDICAL TECH CO., LIMITED



+852 9070 6664



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong