


Mindray DC-85 New Original 3D/4D Curved Array Ultrasound Probe D8-2E

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: Mindray
- Model Number: D8-2E
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification


- Warranty: 90 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: **Mindray DC-85 Ultrasound Probe, 3D/4D Curved Array Ultrasound Probe, D8-2E Ultrasound Probe**

for more products please visit us on 3dultrasoundprobe.com

Product Description

Mindray DC-85 new original 3D/4D curved array ultrasound probe D8-2E

1. Model: D8-2E
2. Application: Obstetrics, Gynaecology, Abdominal
3. Type: 3D/4D Curved array
4. Frequency: 2.0-8.0 Mhz
5. Compatible system: Mindray DC-80,Mindray DC-85



More details welcome to contact with us!

Product Name	Ultrasound Probe/Ultrasound Transducer
Probe Model	Mindray D8-2E
Probe Type	3D/4D Convex array
Central Frequency	2.0-8.0MHz
Compatible system	Mindray DC-80,DC-85
Warranty	90 days
Delivery date	Within 1 week after getting payment
MOQ	1 Unit
Condition	original
Material	Metal and Plastic
Application	Obstetrics, Gynaecology, Abdominal
Service	outright sell/repair

Note: All information provided here is for general information only. Please note that not all functions and features may be available for each system or probe. All information in this publication has been compiled and checked with the utmost diligence, nevertheless errors and mistakes cannot be excluded.

some knowledge about Mindray DC-70:

Mindray's DC-70 ultrasound system is one of Mindray's most advanced ultrasound units available on the market. Adored by diagnosticians, respected by engineers and technicians, and simple to explain to patients, this unit is particularly special. There are many a reason why the DC-70 is so revered; imaging enhancing technology, a multitude of compatible probes / transducers (including single crystal and ComboWave probes) and specialized OB-GYN and women's health applications.

Some of the technology that Mindray has incorporated into this system is the most advanced that they have available. This includes X-engine (to speed up image processing up to 3X the average time), iLive with Hyaline (which provides more clear 3D/4D (live 3D) images than ever before), Smart Planes CNS (enhanced imaging for fetal Central Nervous System), and more.

The imaging technology behind the new available ultrasound probes is staggering. The DC-70 is compatible with single crystal probes - special transducers that implement 3T (Triple-matching layers, Total-cut design, Thermal control) technology to enhance OB-GYN and cardiac ultrasound imaging.. ComboWave transducers, on the other hand, use a new special material to enhance acoustic quality and prevent impediments that occur through echoes. This results in higher quality imaging for small parts (thyroid, breast) and Vascular applications.

Many of the ultrasound transducers that are compatible with many of the probes that also work with the Mindray DC-8 ultrasound system. This makes it easy to rotate between ultrasound units without having to purchase new expensive transducers. The process is easy and fluid.

Ultrasound Features and Specificities:

- Single Crystal Technology Probes
- ComboWaves Technology Probes
- X-Engine
- iLive with Hyaline
- Smart Planes CNS
- Smart Face (fetal enhancing technology)
- Smart FLC
- Smart OB
- Smart NC
- 21.5 inch HD monitor
- 13 inch Touch Screen



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong