

Samsung Medison Cardiac Original Ultrasound System Sonoace X8

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com

Basic Information

- Place of Origin: South Africa
- Brand Name: Samsung Medison
- Model Number: Sonoace X8
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 5-7 days
- Payment Terms: T/T, Western Union
- Supply Ability: 5-10pcs/month



Product Specification

- Warranty: 30 Days
- Lead Time: 5-7 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: **Samsung Medison Cardiac Ultrasound System, Sonoace X8 Ultrasound System**




More Images



Product Description

Samsung Medison cardiac original ultrasound system Sonoace X8



Samsung SonoaceX8 Ultrasound Machine Features

The Samsung Medison SonoAce X8 is a mid-range cart-based ultrasound machine for OB/GYN and cardiac applications. The SonoAce X8 was initially manufactured by Medison starting in 2007 until Samsung acquired the company in 2010.

The SonoACE X8 has several OB/GYN features, including 4D imaging. With 4D, sonographers and clinicians can show expecting parents their unborn child's face and diagnose some conditions before birth, such as cleft lip.

Besides OB/GYN applications, the SonoAce X8 supports several modalities, including cardiac, vascular, MSK, abdominal, and anesthesiology. It has three probe ports and an HD 17-inch LCD monitor, which shows ultrasound scans in high-quality and excellent contrast.

Imaging Modes

2D/3D/4D

Applications

OB/GYN, Vascular, Small Parts, MSK, Anesthesiology, Pediatrics,
Urology, Abdominal

Features

17-inch display, 2D/3D/4D, M-color Flow Mode, Anatomical M-Mode,
Live 3D, Freehand 3D, Color Doppler, Power Angio Doppler, Pulse
Wave Doppler, Trapezoidal Mode, Stress Echo, Strain Rate, Tissue
Harmonic Imaging, Spatial Compounding, Speckle Reduction,
Triplex Mode, Auto IMT



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, KowLoon, Hong Kong