

Mindray Vetus 5 Expert Original Phased Array Ultrasound Probe P8-2

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com

Basic Information

- Place of Origin: China
- Brand Name: Mindray
- Model Number: P8-2
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 90 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: **original phased array ultrasound probe, Mindray Vetus 5 expert Ultrasound Probe, P8-2 Ultrasound Probe**



Product Description

Mindray Vetus5expert original phased array ultrasound probe P8-2

1. Model: P8-2
2. Application: Pediatric, Cardiac, Transcranial, Thoracic, Abdomen,Vascular
3. Type: Phased array
4. Frequency: 2.3-7.2 Mhz
5. Compatible system: Vetus 5 Expert, Vetus 50



More details welcome to contact with us!

Product Name	Ultrasound Probe/Ultrasound Transducer
Probe Model	P8-2
Probe Type	Phased
Central Frequency	2.3-7.2MHz
Compatible system	Vetus 5 Expert, Vetus 50
Warranty	90 days
Delivery date	Within 1 week after getting payment
MOQ	1 Unit
Condition	original
Material	Metal and Plastic
Application	Pediatric, Cardiac, Transcranial, Thoracic, Abdomen,Vascular
Service	Used/new original sell

Mindray Vetus5 Expert

Vetus 5Exp is a cost-effective veterinary ultrasound imaging system with multiple intelligent functions. Configured with brand new VETUS platform, making it easily meet the increasingly diversified veterinarian demands. Incorporating the dedicated animals image presets, simple workflow, and veterinary protective design promise an ideal ultrasound system for your everyday practice.

Key features

Professional ultrasound image presets according to different animals' body characteristics

Animal species are subdivided to different categories by weight and body size for precise diagnosis

Veterinary recognized measurement formula presets, easily user-defined modification

Transducers and the corresponding exam modes are displayed in the same interface, select or switch within 1 second

Recent 4 scanning modes are presented at hand, one-key to switch the frequently-used exam mode

Super smart touch screen supports powerful multi-touch operation
Auto measurement tools – Auto PW, Auto CW maximize the productivity for easy-accurate-scanning
Dedicated application packages – Abdomen/Cardiology/Reproduction/Small Parts
Resourceful tools – vet formula/comment library/animal body mark/report template
iWorks – automated workflow guide iWorks standardize and simplify the vet Abdomen/Cardiology workflow, allowing more focus on the animal diagnosis

Key ergonomic and other features

21.5'' high resolution LED Monitor 178° viewing angle, ultra slim & easy cleaning
13.3'' anti-glare touch screen 30° rotation, easy cleaning
Sliding keypad with protective silicon film, anti-liquid and anti-animal hairs
4 active transducer ports with anti-hairs & dust protective covers
Large capacity battery for long scanning time
300 Mbps Wi-Fi with easy connection everywhere
Anti-dirty multi-functional control panel, rotatable & height adjustable
Unique design for cable management that keep transducer cables off the ground
Stable wheels with small footprint, anti hair-wrapped and 360° rotation



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong