

Samsung Medison Sonoace5500 Curved Array Ultrasound Probe C2-4/30R

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com

Basic Information

- Place of Origin: South Korea
- Brand Name: Samsung Medison
- Model Number: C2-4/30R
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: Samsung Medison Curved array ultrasound probe, Sonoace 5500 Curved array ultrasound probe, C2-4/30R Curved array ultrasound probe



More Images



Product Description

Samsung Medison Sonoace5500 Curved array ultrasound probe C2-4/30R

1. Model: C2-4/30R
2. Application: Pediatrics, Abdominal, Cardiac, Veterinary
3. Type: Curved array
4. Frequency: 2.0-4.0 Mhz
5. Compatible system: Medison SonoAce SA-5500; Medison SonoAce SA-6000C



More details welcome to contact with us!

Product Name	Ultrasound Probe/Ultrasound Transducer
Probe Model	C2-4/30R
Probe Type	Curved array
Central Frequency	2.0-4.0 MHz
Compatible system	Medison SonoAce SA-5500; Medison SonoAce SA-6000C
Warranty	60 days
Delivery date	Within 1 week after getting payment
MOQ	1 Unit
Condition	original
Material	Metal and Plastic
Application	Pediatrics, Abdominal, Cardiac, Veterinary
Service	Used original sale/repair service

Mindray Vetus5 Expert

Vetus 5Exp is a cost-effective veterinary ultrasound imaging system with multiple intelligent functions. Configured with brand new VETUS platform, making it easily meet the increasingly diversified veterinarian demands. Incorporating the dedicated animals image presets, simple workflow, and veterinary protective design promise an ideal ultrasound system for your everyday practice.

Key features

Professional ultrasound image presets according to different animals' body characteristics
 Animal species are subdivided to different categories by weight and body size for precise diagnosis
 Veterinary recognized measurement formula presets, easily user-defined modification
 Transducers and the corresponding exam modes are displayed in the same interface, select or switch within 1 second
 Recent 4 scanning modes are presented at hand, one-key to switch the frequently-used exam mode
 Super smart touch screen supports powerful multi-touch operation
 Auto measurement tools – Auto PW, Auto CW maximize the productivity for easy-accurate-scanning

Dedicated application packages – Abdomen/Cardiology/Reproduction/Small Parts
Resourceful tools – vet formula/comment library/animal body mark/report template
iWorks – automated workflow guide iWorks standardize and simplify the vet Abdomen/Cardiology workflow, allowing more focus on the animal diagnosis

Key ergonomic and other features

21.5" high resolution LED Monitor 178° viewing angle, ultra slim & easy cleaning

13.3" anti-glare touch screen 30° rotation, easy cleaning

Sliding keypad with protective silicon film, anti-liquid and anti-animal hairs

4 active transducer ports with anti-hairs & dust protective covers

Large capacity battery for long scanning time

300 Mbps Wi-Fi with easy connection everywhere

Anti-dirty multi-functional control panel, rotatable & height adjustable

Unique design for cable management that keep transducer cables off the ground

Stable wheels with small footprint, anti hair-wrapped and 360° rotation



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong