

## Canon Toshiba Aplio 500 Endocavity Array Original Ultrasound Probe PVT-781VTE

### Our Product Introduction

#### Basic Information

- Place of Origin: Japan
- Brand Name: Canon Toshiba
- Model Number: PVT-781VTE
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



#### Product Specification

- Warranty: 90 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: **PVT-781VTE Ultrasound Probe, Endocavity Array Ultrasound Probe, Canon Toshiba Aplio 500 Ultrasound Probe**



## Product Description

### Canon Toshiba PVT-781VTE Intravacuity Probe for Aplio 300 Aplio 400 Endocavity Scanner

The Toshiba PVT-781VTE Intra Cavity transducer has a frequency range of 3.6 MHz. - 10.5 MHz. and supports the following applications: Obstetrics, Gynaecology. The Toshiba PVT-781VTE transducer is compatible with: Toshiba Aplio 300 (TUS-A300); Toshiba Aplio 400 (TUS-A400); Toshiba Aplio 500 (TUS-A500) ultrasound systems.

Form Factor: Intra Cavity

Applications: Obstetrics, Gynaecology

Frequency Range: 3.6 - 10.5 MHz.

Compatibility: Toshiba Aplio 300 (TUS-A300); Toshiba Aplio 400 (TUS-A400); Toshiba Aplio 500 (TUS-A500)

Bandwidth: 3.0 – 11.0 MHz

Applications: Transvaginal uterus, ovary

Compatibility: Aplio 300, Aplio 400, Aplio 500

Field of view: 180°

Curvature: Approx. 10 mm

Reusable Biopsy Adapter: UAGV-035A

Fusion Adapter: UA F S -006A





More details welcome to contact with us!

Product Name	Ultrasound Probe/Ultrasound Transducer
Probe Model	Toshiba PVT-781VTE
Probe Type	Endocavity/Transvaginal
Central Frequency	6.5MHz
Compatible system	Toshiba Aplio 300 (TUS-A300)/400 (TUS-A400)/500 (TUS-A500)
Warranty	2-6 months
Delivery date	Within 1 week after getting payment
MOQ	1 Unit
Condition	New or used
Material	Metal and Plastic
Application	Obstetrical, gynecological (OB/GYN) and urological clinical applications and exams.
Elements	128

each system or probe. All information in this publication has been compiled and checked with the utmost diligence, nevertheless errors and mistakes cannot be excluded.

**FOQ:**

Q: What is the MOQ?

A: 1pc.

Q: How's the stock condition?

A: Wide variety of Ultrasound probe/board, new&used&compatible new for your choose.

Q: Is there any guarantee for the goods?

A: Fully-tested and safe packing before shipping, less than 3% warranty return claims of previous orders.

Q: Will you provide warranty?

A: Items are sold with NO WARRANTY, unless it is stated otherwise in the body of the proforma invoice/purchase order. Mostly 90 days for new original, 30/60 days for used.

Q: Do you test all your goods before delivery?

A: Yes, we have 100% tested our goods before delivery.

Q: Can I have the goods returned if it is not working good?

A: Yes, goods are permitted to be returned in same received condition within 7 days. Refund will be handled within 2-4 working days after final confirmation (Freight & bank fee excluded). Shipping fee for return should be borne by buyers

Q: If faults occurred after warranty date, can I contact you?

A: Sure, feel free to contact even years later. We'll look into the faults and try to find out some solutions/suggestions for you reference.

Q: How about the payment terms?

A: All orders are 100% prepaid before delivery. We accept wire transfer/ Western Union from bank, and cash. (US dollars, Euro, or CNY)

Q: How about delivery?

A: We ship worldwide. DHL/EMS/FedEx Express will be arranged within 1-5 working days after payment confirmed unless otherwise specified. We Seller are responsible for export Customs clearance in our end, and buyers for import Customs clearance in their end.



**VOLKA MEDICAL TECH CO., LIMITED**



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong