

Chison IVIS30 IVIS40 IVIS50 IVIS60 Convex Ultrasound Probe D3C60L

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CHISON
- Model Number: D3C60L
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: IVIS50 convex ultrasound probe,
IVIS60 convex ultrasound probe,
IVIS40 convex ultrasound probe



More Images



Product Description

Chison IVIS30 IVIS40 IVIS50 IVIS60 convex ultrasound probe D3C60L

1. Model: D3C60L
2. Application: OB/GYN, Abdominal
3. Type: Convex array
4. Frequency: 3.5MHz
5. Compatible system: Chison IVIS30, Chison IVIS40, Chison IVIS50, Chison IVIS60



More details welcome to contact with us!

Product Name	Chison D3C60L
Probe Type	Convex Array
Frequency	3.5 MHz
Material	Metal and Plastic

Compatible system	Chison IVIS30,Chison IVIS40,Chison IVIS50,Chison IVIS60
Application	OB/GYN,Abdominal
/	/
Condition	Used original,fully tested before delivery.
Warranty	60 days
MOQ	1 unit
available	new&used original,compatible new one

FOQ:

Q: What is the MOQ?

A: 1pc.

Q: How's the stock condition?

A: Wide variety of Ultrasound probe/board, new&used&compatible new for your choose.

Q: Is there any guarantee for the goods?

A: Fully-tested and safe packing before shipping, less than 3% warranty return claims of previous orders.

Q: Will you provide warranty?

A: Items are sold with NO WARRANTY, unless it is stated otherwise in the body of the proforma invoice/purchase order. Mostly 90 days for new original, 30/60 days for used.

Q: Do you test all your goods before delivery?

A: Yes, we have 100% tested our goods before delivery.

Q: Can I have the goods returned if it is not working good?

A: Yes, goods are permitted to be returned in same received condition within 7 days. Refund will be handled within 2-4 working days after final confirmation (Freight & bank fee excluded). Shipping fee for return should be borne by buyers

Q: If faults occurred after warranty date, can I contact you?

A: Sure, feel free to contact even years later. We'll look into the faults and try to find out some solutions/suggestions for you reference.

Q: How about the payment terms?

A: All orders are 100% prepaid before delivery. We accept wire transfer/ Western Union from bank, and cash. (US dollars, Euro, or CNY)

Q: How about delivery?

A: We ship worldwide. DHL/EMS/FedEx Express will be arranged within 1-5 working days after payment confirmed unless otherwise specified. We Seller are responsible for export Customs clearance in our end, and buyers for import Customs clearance in their end.



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong