

ATL IU22 HD11 3D/4D Curved Array Ultrasound Probe V6-2

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com

Basic Information

- Place of Origin: Netherlands
- Brand Name: ATL
- Model Number: V6-2
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: 3D Curved array ultrasound probe,
HD11 Curved array ultrasound probe,
4D Curved array ultrasound probe



More Images



Product Description

ATL IU22 HD11 3D/4D Curved array ultrasound probe V6-2

1. Model: V6-2
2. Application: Abdominal, Obstetrics, Gynaecology, Pediatrics
3. Type: 3D/4D Curved array
4. Frequency: 1.9-4.7 Mhz
5. Compatible system: HD11; HD11XE; IU22; HD5



More details welcome to contact with us!

Tips for ATL IU22 Ultrasound system

Product Features:

- Tissue Harmonic Imaging with Pulse Inversion technology
- Harmonic SonoCT imaging
- XRES adaptive image processing technology
- iSCAN intelligent scanning for one-button TGC, gain and compression map optimization
- iSCAN with AGC for real-time line-by-line TGC optimization
- Simultaneous 2D M-mode
- Broadband Flow imaging
- CPA
- 2D/PW Doppler
- CW/PW
- Contrast Specific Imaging (CSI)
- Panoramic imaging with SonoCT, XRES and harmonic modes
- Freehand 3D and MPR imaging with SonoCT, XRES and harmonic modes
- Automated 3D and 4D imaging with SonoCT, XRES and harmonic modes
- Chroma imaging in 2D, 3D, 4D, MPR, Panoramic, M-mode and Doppler modes
- DICOM Print/Store
- Ergonomics human-center design for comfort and convenience

ATL IU22 Applications:

Abdominal

Obstetrical

Vascular

Urology

Cerebrovascular

Peripheral vascular

Gynecological and fertility

Small parts and superficial

Musculoskeletal

Pediatrics general imaging

Prostate

Adult cardiology

Contrast Imaging



VOLKA MEDICAL TECH CO., LIMITED



+8619536805795



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong