

Repair Philips C6-2 Affiniti 50 Convex OB/GYN Probe C6-2

Our Product Introduction

Basic Information

- Place of Origin: USA
- Brand Name: Philips
- Model Number: C6-2
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 working days
- Payment Terms: T/T, Western Union
- Supply Ability: 30-50pcs/month



Product Specification

- Warranty: 60 Days After Fixed
- Delivery Method: Express Or As Clients Required
- Before Shipping: Test Images&videos Shared
- Service: Repair
- Highlight: **OB/GYN Probe Repair , C6-2 Probe Repair**

Product Description

Repair Philips C6-2 Affiniti 50 Convex OB/GYN Probe C6-2

1. Model: C6-2
2. Type: Convex
3. Application: OB/GYN, Abdominal
4. Compatible system: Affiniti 50, CX30, Epiq, Sparq
5. Repair including: replace crystal, repair motor control problem, replace dome and etc
6. With 60 days warranty after fixed and 3-5 working days to repair.



How to do complete clinical evaluation of your transducer?

1. Frequency conversion

Whether the frequency can be transformed, it's the important indicator.

2. Frequency range

The wider the range, the better.

3. Scanning angle

The bigger the angle, the better. Especially for the intracavitary probes.

