

Canon Toshiba Aplio 300 Sector Cardiac Transcranial Ultrasound Probe PST-25BT

Our Product Introduction

for more products please visit us on 3dultrasoundprobe.com


Basic Information

- Place of Origin: Japan
- Brand Name: Canon Toshiba
- Model Number: PST-25BT
- Minimum Order Quantity: 1pc
- Price: negotiable
- Packaging Details: carton
- Delivery Time: 3-5 days
- Payment Terms: T/T, Western Union
- Supply Ability: 10-30pcs/month



Product Specification

- Warranty: 60 Days
- Lead Time: 3-5 Days
- Service: Outright
- Shipping Method: Express, or As Clients Required
- Highlight: Canon Toshiba Cardiac Transcranial ultrasound probe
,
PST-25BT Cardiac Transcranial ultrasound probe



Product Description

Canon Toshiba Aplio 300 Sector Cardiac Transcranial Ultrasound Probe PST-25BT

Product Specifications

Attribute	Value
Warranty	60 days
Lead time	3-5 days
Service	Outright
Shipping method	Express, or as clients required



Detailed Product Information

Model: PST-25BT

Application: Cardiac, Transcranial

Type: Sector

Frequency: 1.8/4.2 MHz

Compatible Systems: Aplio 300, Aplio 50 SSA-750A, Aplio 500, Aplio MX SSA-780A, Aplio XG SSA-790A, Xario SSA-660A

Delivery: Within 1 week after payment

MOQ: 1 Unit

Condition: Original

Material: Metal and Plastic

Clinical Performance Features

The Canon Toshiba Aplio 300 PST-25BT ultrasound probe delivers outstanding image quality for cardiac and transcranial applications, supporting confident diagnoses with:

Live iSlice technology for precise region of interest examination

Advanced tissue characterization capabilities

Optimized imaging for specific exam types

Reliable performance in clinical settings

Ultrasound Applications

Abdominal

Cardiac

Vascular

Small Parts

Transcranial

OB-GYN

Pediatrics

Compatible Ultrasound Probes and Transducers

ATL S4-2 sector array probe

ATL L12-5 linear array probe

ATL L9-3 linear array probe

ATL C6-3 curved array probe

ATL C5-2 curved array probe

ATL C8-5 curved array probe

ATL V9-4V endocavity probe

For more details or inquiries, please contact us.



VOLKA MEDICAL TECH CO., LIMITED



+852 9070 6664



info@volka-medtech.com



3dultrasoundprobe.com

Office 5, 8/F, Mega Cube, 8 Wang Kwong Road, Kowloon Bay, Kowloon, Hong Kong